

Product datasheet for **HP218141**

LYSMD2 Human qPCR Primer Pair (NM_153374)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 256586 |
| Forward Sequence: | CTGAGGAAGTGTCAGCCAGAGA |
| Reverse Sequence: | GTGATAGAGGGAAGTTGCATAGG |
| ACCN: | BC022075 , NM_153374 , NM_153374.1 , NM_153374.2 , BC022075.1 , BC033515 , BC033515.1 , BI822634 , BQ446302 , NM_153374.3 |
| UniProt ID: | Q8IV50 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)