

## Product datasheet for **HP218146**

### Killer cell immunoglobulin like receptor 3DL3 (KIR3DL3) Human qPCR Primer Pair (NM\_153443)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	115653
Forward Sequence:	TCCTTCATCGCTGGTGTGCCAA
Reverse Sequence:	CAACTGTGCGTATGTACCTCC
ACCN:	<u>NM_153443</u> , <u>NM_153443.1</u> , <u>NM_153443.2</u> , <u>NM_153443.3</u> , <u>NM_153443.4</u> , <u>BC136283</u> , <u>BC143802</u>
Synonyms:	CD158Z; KIR3DL7; KIR44; KIRC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)