

Product datasheet for **HP218166**

ZNF114 Human qPCR Primer Pair (NM_153608)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	163071
Forward Sequence:	CTCTACAGAGACGTGATGCTGG
Reverse Sequence:	CACACTCTGTTGGCTTCAGGAAA
ACCN:	<u>NM_153608</u> , <u>NM_153608.1</u> , <u>NM_153608.2</u> , <u>NM_153608.3</u> , <u>BC125069</u> , <u>BC014935</u> , <u>BC016811</u> , <u>BC034032</u> , <u>BC050996</u> , <u>BC125068</u> , <u>BG476217</u> , <u>BU182389</u> , <u>NM_153608.4</u>
UniProt ID:	<u>Q8NC26</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)