

Product datasheet for **HP218168**

CYB561A3 Human qPCR Primer Pair (NM_153611)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	220002
Forward Sequence:	CTCTGTCCATCGCATCCGTCAT
Reverse Sequence:	TGCTGTTGGCAAAGACCGCCTC
ACCN:	NM_153611 , NM_153611.1 , NM_153611.2 , NM_153611.3 , NM_153611.4 , BC004391 , BC004391.1 , BC014045 , BC047710 , BM692535 , BP296553 , BP366144 , BP396674
UniProt ID:	Q8NBI2
Synonyms:	CYBASC3; LCYTB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)