

## Product datasheet for **HP218170**

### AGPAT7 (LPCAT4) Human qPCR Primer Pair (NM\_153613)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	254531
Forward Sequence:	GCTACAGCTCTCTGATCCTCAG
Reverse Sequence:	GACGAGTTAGCTCTCCAGGCT
ACCN:	<a href="#">NM_153613</a> , <a href="#">NM_153613.1</a> , <a href="#">NM_153613.2</a> , <a href="#">BC092463</a> , <a href="#">BC024892</a> , <a href="#">BC028209</a> , <a href="#">BC048430</a> , <a href="#">BC144236</a> , <a href="#">NM_153613.3</a>
UniProt ID:	<a href="#">Q643R3</a>
Synonyms:	AGPAT7; AYTL3; LPAAT-eta; LPEAT2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)