

Product datasheet for **HP218178**

ZNF475 (ZFP1) Human qPCR Primer Pair (NM_153688)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	162239
Forward Sequence:	CCTATACATGGATGTGATGCTGG
Reverse Sequence:	TCGTCTGGGTTTCTGTGGTCTC
ACCN:	NM_153688 , NM_153688.1 , NM_153688.2 , NM_153688.3 , BC033774 , BC064604 , BC113863 , BC114471 , BC127184 , BC160028 , BP219674 , BX647144 , NM_153688.4
UniProt ID:	Q6P2D0
Synonyms:	ZNF475
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)