

Product datasheet for **HP218212**

C10orf46 (CACUL1) Human qPCR Primer Pair (NM_153810)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	143384
Forward Sequence:	GCCCAGTCAACACCATTTCAGG
Reverse Sequence:	AGAGTTGGAGCCATCTGAACCC
ACCN:	<u>NM_153810</u> , <u>NM_153810.1</u> , <u>NM_153810.2</u> , <u>NM_153810.3</u> , <u>NM_153810.4</u> , <u>BC017894</u> , <u>BC038294</u> , <u>BG772719</u> , <u>BU623607</u> , <u>NM_153810.5</u>
UniProt ID:	<u>Q86Y37</u>
Synonyms:	C10orf46; CAC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)