

Product datasheet for **HP218259**

SALF (STON1-GTF2A1L) Human qPCR Primer Pair (NM_172311)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	286749
Forward Sequence:	TATGAGAGTGCCTACCAGGCAG
Reverse Sequence:	ACTGAACAGTAGCAAATGGATACC
ACCN:	<u>NM_172311</u> , <u>NM_172311.1</u> , <u>NM_172311.2</u> , <u>BC101772</u> , <u>BC101772.1</u> , <u>BC143520</u> , <u>BC101774</u> , <u>BG191372</u> , <u>BX091338</u> , <u>NM_172311.3</u>
UniProt ID:	<u>Q9Y6Q2</u>
Synonyms:	ALF; GTF2A1L; GTF2A1LF; SALF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)