

## Product datasheet for **HP218274**

### ARMC3 Human qPCR Primer Pair (NM\_173081)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	219681
Forward Sequence:	AACTGCGTGTGACGTTGAAGCC
Reverse Sequence:	AGCACAGACCATCACTGCCCAA
ACCN:	<u><a href="#">NM_173081</a></u> , <u><a href="#">NM_173081.1</a></u> , <u><a href="#">NM_173081.2</a></u> , <u><a href="#">NM_173081.3</a></u> , <u><a href="#">NM_173081.4</a></u> , <u><a href="#">BC039312</a></u> , <u><a href="#">BC039312.1</a></u> , <u><a href="#">BX115364</a></u>
UniProt ID:	<u><a href="#">Q5W041</a></u>
Synonyms:	CT81; KU-CT-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)