

Product datasheet for **HP218275**

LIN9 Human qPCR Primer Pair (NM_173083)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	286826
Forward Sequence:	GCCTTTGGACAAAAACAGCGGC
Reverse Sequence:	TTCTCCACGGCGACTGTCCTAA
ACCN:	NM_173083 , NM_173083.1 , NM_173083.2 , NM_173083.3 , BC045625 , BC045625.1 , BC043444 , BC065302 , BX537869
UniProt ID:	Q5TKA1
Synonyms:	BARA; BARPsv; Lin-9; TGS; TGS1; TGS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)