

Product datasheet for **HP218277**

Cytokeratin 6C (KRT6C) Human qPCR Primer Pair (NM_173086)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	286887
Forward Sequence:	AAGCCAAGGCAGACTCTCAC
Reverse Sequence:	TGGATAGCACCCGGATGTGTC
ACCN:	<u>NM_173086</u> , <u>NM_173086.1</u> , <u>NM_173086.2</u> , <u>NM_173086.3</u> , <u>NM_173086.4</u> , <u>BC110639</u> , <u>BC130583</u> , <u>BC130585</u>
UniProt ID:	<u>P48668</u>
Synonyms:	K6E; KRT6E; PPKNEFD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)