

Product datasheet for **HP218294**

MOB4A (MOB1B) Human qPCR Primer Pair (NM_173468)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	92597
Forward Sequence:	TTCGGATGGCTGTCATGCTTCC
Reverse Sequence:	GCTGACATCACTGGACAACCTCTC
ACCN:	<u>NM_173468</u> , <u>NM_173468.1</u> , <u>NM_173468.2</u> , <u>NM_173468.3</u> , <u>BC038112</u> , <u>BC038112.1</u> , <u>BC044913</u> , <u>BF129919</u> , <u>BG033976</u> , <u>BG433921</u> , <u>BQ229138</u> , <u>BX089615</u> , <u>NM_173468.4</u>
UniProt ID:	<u>Q7L9L4</u>
Synonyms:	MATS2; MOB4A; MOBKL1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)