

## Product datasheet for **HP218348**

### C18orf54 Human qPCR Primer Pair (NM\_173529)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	162681
Forward Sequence:	ATCTCCTGTTCCCGTTAACTCTG
Reverse Sequence:	CAGGACCATGATGTTTGCCTCC
ACCN:	<a href="#">NM_173529</a> , <a href="#">NM_173529.1</a> , <a href="#">NM_173529.2</a> , <a href="#">NM_173529.3</a> , <a href="#">NM_173529.4</a> , <a href="#">NM_173529.5</a> , <a href="#">BC036054</a> , <a href="#">BC036054.2</a> , <a href="#">BC043651</a> , <a href="#">BX640877</a> , <a href="#">NM_173529.6</a>
UniProt ID:	<a href="#">Q8IYD9</a>
Synonyms:	LAS2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)