

Product datasheet for **HP218374**

RFXDC1 (RFX6) Human qPCR Primer Pair (NM_173560)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222546
Forward Sequence:	ACTCACTGCCAGTGTATCCTGG
Reverse Sequence:	GGAAGGAGATGGTCAGGCATTC
ACCN:	NM_173560 , NM_173560.1 , NM_173560.2 , NM_173560.3 , BC039248 , BC039248.1 , BQ269808
UniProt ID:	Q8HWS3
Synonyms:	dj955L16.1; MTCHRS; MTF5; RFXDC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)