

Product datasheet for **HP218414**

INO80E Human qPCR Primer Pair (NM_173618)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	283899
Forward Sequence:	CTCAAGTTCCTCATCTACGAGCA
Reverse Sequence:	GTCCGAAGAGTCTTCATCCACG
ACCN:	NM_173618 , NM_173618.1 , NM_173618.2 , BC047712 , BC047712.2 , BC035693 , BI093312 , BM924992 , BU528528 , BX356644 , NM_173618.3
UniProt ID:	Q8NBZ0
Synonyms:	CCDC95
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)