

Product datasheet for **HP218510**

CES4A Human qPCR Primer Pair (NM_173815)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	283848
Forward Sequence:	GCGAGAAAGTGGTCTGGTGTT
Reverse Sequence:	AAGGCTGCGATGTTCTCTCTGCA
ACCN:	<u>NM_173815</u> , <u>NM_173815.2</u> , <u>NM_173815.3</u> , <u>NM_173815.4</u> , <u>NM_173815.5</u> , <u>NM_173815.6</u> , <u>BC029039</u> , <u>BC041823</u> , <u>BC064573</u> , <u>BC084555</u> , <u>BC166638</u> , <u>BI914453</u> , <u>NM_173815.7</u>
UniProt ID:	<u>Q5XG92</u>
Synonyms:	CES6; CES8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)