

## Product datasheet for **HP218518**

### SFRS12IP1 (SREK1IP1) Human qPCR Primer Pair (NM\_173829)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	285672
Forward Sequence:	AATTGCAGGCATTACAGGAAAAAAG
Reverse Sequence:	GTGTCCTCTTCAGTGGAACTGG
ACCN:	<a href="#">NM_173829</a> , <a href="#">NM_173829.1</a> , <a href="#">NM_173829.2</a> , <a href="#">NM_173829.3</a> , <a href="#">BC127085</a> , <a href="#">BC047006</a> , <a href="#">BC108716</a> , <a href="#">BC127086</a> , <a href="#">NM_173829.4</a>
UniProt ID:	<a href="#">Q8N9Q2</a>
Synonyms:	P18SRP; SFRS12IP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)