

Product datasheet for **HP218538**

ANUBL1 (ZFAND4) Human qPCR Primer Pair (NM_174890)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	93550
Forward Sequence:	GGTGTACCTTGAAGCTAGTTTTGG
Reverse Sequence:	GACCTCAACTCGACTGGAATCC
ACCN:	<u>NM_174890</u> , <u>NM_174890.1</u> , <u>NM_174890.2</u> , <u>NM_174890.3</u> , <u>BC045587</u> , <u>BC017591</u> , <u>BC048968</u> , <u>NM_174890.4</u>
UniProt ID:	<u>Q86XD8</u>
Synonyms:	ANUBL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)