

Product datasheet for **HP218548**

FAM9C Human qPCR Primer Pair (NM_174901)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	171484
Forward Sequence:	GGAAAGAAAACCTGTTACAGAGAC
Reverse Sequence:	GCTGCAATATCTTCTAGCTCCTTG
ACCN:	<u>NM_174901</u> , <u>NM_174901.1</u> , <u>NM_174901.2</u> , <u>NM_174901.3</u> , <u>NM_174901.4</u> , <u>NM_174901.5</u> , <u>BC127957</u> , <u>BX647580</u> , <u>BX648111</u> , <u>BX648919</u> , <u>NM_174901.6</u>
UniProt ID:	<u>Q8IZT9</u>
Synonyms:	TEX39C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)