

Product datasheet for **HP218551**

FAM98C Human qPCR Primer Pair (NM_174905)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	147965
Forward Sequence:	CATCTGCTTTCCACTGGAGTGAC
Reverse Sequence:	GGCAGCCAGAACGTGTGCAATG
ACCN:	<u>BC034383</u> , <u>NM_174905</u> , <u>NM_174905.1</u> , <u>NM_174905.2</u> , <u>NM_174905.3</u> , <u>BC117258</u> , <u>BC036482</u> , <u>BC080606</u> , <u>BC100033</u> , <u>BC143830</u> , <u>BC171747</u> , <u>BI871449</u> , <u>BM561971</u> , <u>BX641942</u>
UniProt ID:	<u>Q17RN3</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)