

Product datasheet for **HP218563**

SH3D20 (ARHGAP27) Human qPCR Primer Pair (NM_174919)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	201176
Forward Sequence:	AAGCAGCGGCAGCTTCAAGGC
Reverse Sequence:	CAAGTCCTGGATGACCTCGTAG
ACCN:	NM_174919 , NM_174919.3 , BC067345 , BC101388 , BC101391 , BC029634 , BC101389 , BC101390 , BX386761
UniProt ID:	Q6ZUM4
Synonyms:	CAMGAP1; PP905; SH3D20; SH3P20
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)