

Product datasheet for **HP218566**

ADCK5 Human qPCR Primer Pair (NM_174922)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	203054
Forward Sequence:	GGCTCTACCAGTTCCTGGAGG
Reverse Sequence:	ATCTCGGCGAACAGGAGGTAGT
ACCN:	<u>NM_174922</u> , <u>NM_174922.1</u> , <u>NM_174922.2</u> , <u>NM_174922.3</u> , <u>NM_174922.4</u> , <u>BC101659</u> , <u>BC031570</u> , <u>BC032402</u> , <u>BC064345</u> , <u>BC085013</u> , <u>BC143322</u> , <u>BC143326</u> , <u>BI825038</u> , <u>BM146288</u> , <u>BU845476</u> , <u>NM_174922.5</u>
UniProt ID:	<u>Q3MIX3</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)