

Product datasheet for **HP218581**

CD163L1 Human qPCR Primer Pair (NM_174941)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	283316
Forward Sequence:	GCCGAAAAGTTCCAGTGTGAAGG
Reverse Sequence:	GCCATTCACAAGTCGGACATCTG
ACCN:	<u>NM_174941</u> , <u>NM_174941.1</u> , <u>NM_174941.2</u> , <u>NM_174941.3</u> , <u>NM_174941.4</u> , <u>BC104962</u> , <u>BC104964</u> , <u>BQ024006</u>
UniProt ID:	<u>Q9NR16</u>
Synonyms:	CD163B; M160; SCARI2; WC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)