

Product datasheet for **HP218593**

SEC14L3 Human qPCR Primer Pair (NM_174975)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	266629
Forward Sequence:	CGGCTCATCACACCAAGTGGAA
Reverse Sequence:	GGAAAACCTCCGAAGCCGATGTC
ACCN:	<u>NM_174975</u> , <u>NM_174975.1</u> , <u>NM_174975.2</u> , <u>NM_174975.3</u> , <u>NM_174975.4</u> , <u>BC101004</u> , <u>BC069641</u> , <u>BC101001</u> , <u>BC101002</u> , <u>BC101003</u> , <u>BE219833</u> , <u>BQ186227</u> , <u>BU584962</u> , <u>NM_174975.5</u>
UniProt ID:	<u>Q9UDX4</u>
Synonyms:	TAP2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)