

## Product datasheet for **HP218601**

### H4-16 Human qPCR Primer Pair (NM\_175054)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	121504
Forward Sequence:	CGTCAAGCGCATTCTGGTCTC
Reverse Sequence:	TTTCAGCGGTACACCACATCC
ACCN:	<u><a href="#">BC020884</a></u> , <u><a href="#">NM_175054</a></u> , <u><a href="#">NM_175054.1</a></u> , <u><a href="#">NM_175054.2</a></u> , <u><a href="#">BC111093</a></u>
UniProt ID:	<u><a href="#">P62805</a></u>
Synonyms:	H4/p; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; HIST4H4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)