

Product datasheet for **HP218610**

H2AC21 Human qPCR Primer Pair (NM_175065)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	317772
Forward Sequence:	GGGACAACAAGAAGACGCGCAT
Reverse Sequence:	CAGCCTGGATATTGGCAAGAC
ACCN:	<u>NM_175065</u> , <u>NM_175065.1</u> , <u>NM_175065.2</u> , <u>BC172330</u> , <u>NM_175065.3</u>
UniProt ID:	<u>Q8IUE6</u>
Synonyms:	H2AB; HIST2H2AB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)