

## Product datasheet for **HP218642**

### Alpha Taxilin (TXLNA) Human qPCR Primer Pair (NM\_175852)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	200081
Forward Sequence:	AGCTTGAGAGCCTATGCCGTGA
Reverse Sequence:	GTCATTGTGCTGTTCCATCTGC
ACCN:	<a href="#">NM_175852</a> , <a href="#">NM_175852.1</a> , <a href="#">NM_175852.2</a> , <a href="#">NM_175852.3</a> , <a href="#">BC080578</a> , <a href="#">BC029686</a> , <a href="#">BC046565</a> , <a href="#">BC065263</a> , <a href="#">BC103823</a> , <a href="#">BC103824</a> , <a href="#">BC103825</a> , <a href="#">NM_175852.4</a>
UniProt ID:	<a href="#">P40222</a>
Synonyms:	IL14; TXLN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)