

Product datasheet for **HP218696**

AGR3 Human qPCR Primer Pair (NM_176813)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	155465
Forward Sequence:	GCCTCCTCAGACTCTCAAGA
Reverse Sequence:	GAGAGATTGACAATCCTCCAGG
ACCN:	<u>NM_176813</u> , <u>NM_176813.1</u> , <u>NM_176813.2</u> , <u>NM_176813.4</u> , <u>BC058284</u> , <u>BC058284.1</u> , <u>BG540617</u> , <u>BM984550</u> , <u>BU682778</u> , <u>NM_176813.5</u>
UniProt ID:	<u>Q8TD06</u>
Synonyms:	AG-3; AG3; BCMP11; hAG-3; HAG3; PDIA18
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)