

Product datasheet for **HP218697**

ZNF800 Human qPCR Primer Pair (NM_176814)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	168850
Forward Sequence:	TTCTACTGCCACCAAGTCTCC
Reverse Sequence:	TTCTCGTTTGCCACACTTGATA
ACCN:	NM_176814 , NM_176814.1 , NM_176814.2 , NM_176814.4 , BC110623 , BC037348 , BC047641 , BC110624 , BG112338 , BU581251 , BX537723 , NM_176814.5
UniProt ID:	Q2TB10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)