

Product datasheet for **HP218705**

NICE2 (S100A7A) Human qPCR Primer Pair (NM_176823)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	338324
Forward Sequence:	CTGGTCAGTCTCCAGCTAACCT
Reverse Sequence:	TGTGGAGTGGAGCAGTCATTCC
ACCN:	<u>NM_176823</u> , <u>NM_176823.1</u> , <u>NM_176823.2</u> , <u>NM_176823.3</u> , <u>BC040254</u> , <u>BC148486</u> , <u>BC153082</u> , <u>NM_176823.4</u>
UniProt ID:	<u>Q86SG5</u>
Synonyms:	NICE-2; NICE2; S100A7f; S100A7L1; S100A15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)