

Product datasheet for **HP218735**

HCAR2 Human qPCR Primer Pair (NM_177551)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	338442
Forward Sequence:	CTGTTTCCACCTCAAGTCCTGG
Reverse Sequence:	GCCTCACATAGTTGTCCATCAGG
ACCN:	<u>NM_177551</u> , <u>NM_177551.1</u> , <u>NM_177551.2</u> , <u>NM_177551.3</u> , <u>BC027965</u> , <u>BC027965.1</u> , <u>BC056419</u> , <u>BC063461</u> , <u>NM_177551.4</u>
UniProt ID:	<u>Q8TDS4</u>
Synonyms:	GPR109A; HCA2; HM74a; HM74b; NIACR1; Puma-g; PUMAG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)