

## Product datasheet for **HP218739**

### PDE12 Human qPCR Primer Pair (NM\_177966)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	201626
Forward Sequence:	AGGAGAAGGTGCTCCAGAGATC
Reverse Sequence:	GGCTACTGCCATTTGAATGAGGC
ACCN:	<a href="#">NM_177966</a> , <a href="#">NM_177966.1</a> , <a href="#">NM_177966.2</a> , <a href="#">NM_177966.3</a> , <a href="#">NM_177966.4</a> , <a href="#">NM_177966.5</a> , <a href="#">NM_177966.6</a> , <a href="#">BC034978</a> , <a href="#">BM713922</a> , <a href="#">BX647169</a> , <a href="#">NM_177966.7</a>
UniProt ID:	<a href="#">Q6L8Q7</a>
Synonyms:	2'-PDE; 2-PDE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)