

Product datasheet for **HP218752**

C1orf201 (STPG1) Human qPCR Primer Pair (NM_178122)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 90529 |
| Forward Sequence: | CATTCGTGTCCAATACCAGCCG |
| Reverse Sequence: | GTGTGACATCCCTACAGAACCG |
| ACCN: | NM_178122 , NM_178122.1 , NM_178122.2 , NM_178122.4 , BC047705 , BC047705.1 , BC035061 , BC017650 , BC034315 , BC040212 , BC063891 , BI917854 , BM922375 |
| UniProt ID: | Q5TH74 |
| Synonyms: | C1orf201; MAPO2 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)