

## Product datasheet for **HP218761**

### CYP4Z1 Human qPCR Primer Pair (NM\_178134)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	199974
Forward Sequence:	CTGCTTGGCAAAGTACCCTGAG
Reverse Sequence:	GTGGTGTAAGGCATCTGGCTCA
ACCN:	<u><a href="#">NM_178134</a></u> , <u><a href="#">NM_178134.1</a></u> , <u><a href="#">NM_178134.2</a></u> , <u><a href="#">BC140382</a></u> , <u><a href="#">BC146466</a></u> , <u><a href="#">NM_178134.3</a></u>
UniProt ID:	<u><a href="#">Q86W10</a></u>
Synonyms:	CYP4A20
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)