

Product datasheet for **HP218771**

GPIHBP1 Human qPCR Primer Pair (NM_178172)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	338328
Forward Sequence:	AGGTGGAAGAGGAGGAGACCAA
Reverse Sequence:	GTGAGCAGTTCTGCGTCAGGTT
ACCN:	NM_178172 , NM_178172.1 , NM_178172.2 , NM_178172.3 , NM_178172.4 , NM_178172.5 , BC035810 , BC063857
UniProt ID:	Q8IV16
Synonyms:	GPI-HBP1; HYPL1D
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)