

## Product datasheet for **HP218782**

### eEF1A2 binding protein (IGFN1) Human qPCR Primer Pair (NM\_178275)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	91156
Forward Sequence:	CAACACCTGGACGGCAGTGAAC
Reverse Sequence:	GCAAAGACAGCATCCGATGGAG
ACCN:	<u><a href="#">NM_178275</a></u> , <u><a href="#">NM_178275.1</a></u> , <u><a href="#">NM_178275.2</a></u> , <u><a href="#">NM_178275.3</a></u> , <u><a href="#">NM_178275.4</a></u> , <u><a href="#">BC131529</a></u> , <u><a href="#">BC028184</a></u> , <u><a href="#">BC033185</a></u>
UniProt ID:	<u><a href="#">Q86VF2</a></u>
Synonyms:	DKFZp434B1231; eEF1A2 binding protein; EEF1A2BP1; immunoglobulin-like and fibronectin type III domain containing 1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)