

Product datasheet for **HP218785**

RILPL1 Human qPCR Primer Pair (NM_178314)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	353116
Forward Sequence:	GGAGGAGAACAAGCAGCTCATG
Reverse Sequence:	GTTTGTCACCACCTCCTCAG
ACCN:	<u>NM_178314</u> , <u>NM_178314.1</u> , <u>NM_178314.2</u> , <u>NM_178314.3</u> , <u>NM_178314.4</u> , <u>BC041375</u> , <u>BC080626</u> , <u>BC089444</u> , <u>BX281699</u> , <u>BX359997</u> , <u>NM_178314.5</u>
UniProt ID:	<u>Q5EBL4</u>
Synonyms:	GOSPEL; RLP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)