

Product datasheet for **HP218787**

MRPL52 Human qPCR Primer Pair (NM_178336)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	122704
Forward Sequence:	GAGCTCCCAGACTGGTCATATG
Reverse Sequence:	CAGCGTCCATTTCTGTGACAG
ACCN:	NM_178336 , NM_178336.1 , NM_178336.2 , BC068070 , BC068070.1 , BC055381 , BG714882 , BG716120 , BQ082981 , BU600884 , BX248749 , NM_178336.3
UniProt ID:	Q86TS9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)