

## Product datasheet for **HP218808**

### **MARCHF3 Human qPCR Primer Pair (NM\_178450)**

#### **Product data:**

Product Type:	qPCR Primer Pairs
Gene ID:	115123
Forward Sequence:	CTGGCTGTCATCCTCAAACACC
Reverse Sequence:	TGTCGCCAAACAGAGTCCGCTT
ACCN:	<u><a href="#">NM_178450</a></u> , <u><a href="#">NM_178450.1</a></u> , <u><a href="#">NM_178450.2</a></u> , <u><a href="#">NM_178450.3</a></u> , <u><a href="#">NM_178450.4</a></u> , <u><a href="#">BC047569</a></u> , <u><a href="#">BC047569.1</a></u> , <u><a href="#">BC146948</a></u> , <u><a href="#">BC146964</a></u> , <u><a href="#">BG719261</a></u> , <u><a href="#">NM_178450.5</a></u>
UniProt ID:	<u><a href="#">Q86UD3</a></u>
Synonyms:	MARCH-III; MARCH3; RNF173
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)