

## Product datasheet for **HP218809**

### LRR50 (DNAAF1) Human qPCR Primer Pair (NM\_178452)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	123872
Forward Sequence:	CCAGTGTTTCCAAAGGACAGAGC
Reverse Sequence:	CCAAGGCTTCAATGCTGTCTGTG
ACCN:	<a href="#">NM_178452</a> , <a href="#">NM_178452.1</a> , <a href="#">NM_178452.2</a> , <a href="#">NM_178452.3</a> , <a href="#">NM_178452.4</a> , <a href="#">NM_178452.5</a> , <a href="#">BC024009</a> , <a href="#">BC024009.2</a> , <a href="#">BX648657</a> , <a href="#">NM_178452.6</a>
UniProt ID:	<a href="#">Q8NEP3</a>
Synonyms:	CILD13; DAU1; LRR50; ODA7; swt
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)