

## Product datasheet for **HP218905**

### **IQUB Human qPCR Primer Pair (NM\_178827)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	154865
<b>Forward Sequence:</b>	CTCCCAAGATGAGAGGCTGGAT
<b>Reverse Sequence:</b>	CTCCTCTCATCATAAGGTCAACC
<b>ACCN:</b>	<a href="#">NM_178827</a> , <a href="#">NM_178827.1</a> , <a href="#">NM_178827.2</a> , <a href="#">NM_178827.3</a> , <a href="#">NM_178827.4</a> , <a href="#">BC128181</a> , <a href="#">BC026173</a> , <a href="#">BC032838</a> , <a href="#">BC036121</a> , <a href="#">BC057256</a> , <a href="#">BC068607</a> , <a href="#">BC091520</a> , <a href="#">BC128182</a> , <a href="#">NM_178827.5</a>
<b>UniProt ID:</b>	<a href="#">Q8NA54</a>
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)