

Product datasheet for **HP218930**

GPR180 Human qPCR Primer Pair (NM_180989)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	160897
Forward Sequence:	GATTCTACGCCTGCATCCACTG
Reverse Sequence:	GATCCCTGCTAAGTTGTGGTGTG
ACCN:	NM_180989 , NM_180989.1 , NM_180989.2 , NM_180989.3 , NM_180989.4 , NM_180989.5 , BC052243 , BC052243.1 , BC033906 , BG112926 , BI856811 , BX348762
UniProt ID:	Q86V85
Synonyms:	ITR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)