

## Product datasheet for **HP218978**

### PIFO Human qPCR Primer Pair (NM\_181643)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	128344
Forward Sequence:	ACTTTGGGAGCCAGAACAGCTG
Reverse Sequence:	CCTGCTGAGGACTTACTTTCTGG
ACCN:	<a href="#">NM_181643</a> , <a href="#">NM_181643.1</a> , <a href="#">NM_181643.2</a> , <a href="#">NM_181643.3</a> , <a href="#">NM_181643.4</a> , <a href="#">NM_181643.5</a> , <a href="#">BC029237</a> , <a href="#">BC050319</a> , <a href="#">BC068603</a> , <a href="#">BC093891</a> , <a href="#">BC101501</a> , <a href="#">BF674731</a> , <a href="#">BU686257</a> , <a href="#">NM_181643.6</a>
UniProt ID:	<a href="#">Q8TCI5</a>
Synonyms:	C1orf88; pitchfork
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)