

Product datasheet for **HP218986**

Leukotriene B4 Receptor (LTB4R) Human qPCR Primer Pair (NM_181657)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1241
Forward Sequence:	CCTGTGTCACCTATGTCTGCGGA
Reverse Sequence:	ATGCCTTGGTGCCTAGCTTCT
ACCN:	<u>NM_181657</u> , <u>NM_181657.1</u> , <u>NM_181657.2</u> , <u>NM_181657.3</u> , <u>BC004545</u> , <u>BQ182452</u> , <u>BT007267</u> , <u>BU146237</u>
UniProt ID:	<u>Q15722</u>
Synonyms:	BLT1; BLTR; CMKRL1; GPR16; LTB4R1; LTBR1; P2RY7; P2Y7
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)