

## Product datasheet for **HP219012**

### GPR120 (FFAR4) Human qPCR Primer Pair (NM\_181745)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	338557
Forward Sequence:	ACCTCGGAACACCTCCTGGATG
Reverse Sequence:	ATCTGGTGGCTCTCCGAGTAGG
ACCN:	<u><a href="#">NM_181745</a></u> , <u><a href="#">NM_181745.1</a></u> , <u><a href="#">NM_181745.2</a></u> , <u><a href="#">NM_181745.3</a></u> , <u><a href="#">BC101175</a></u> , <u><a href="#">BM757151</a></u> , <u><a href="#">BX647721</a></u> , <u><a href="#">NM_181745.4</a></u>
UniProt ID:	<u><a href="#">Q5NUL3</a></u>
Synonyms:	BMIQ10; GPR120; GPR129; GT01; O3FAR1; PGR4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)