

Product datasheet for **HP219067**

CBLN2 Human qPCR Primer Pair (NM_182511)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	147381
Forward Sequence:	CAACCACTTTGATCTTGCTTCCAG
Reverse Sequence:	CCGAGATCACTGGGTAGCCATT
ACCN:	<u>NM_182511</u> , <u>NM_182511.1</u> , <u>NM_182511.2</u> , <u>NM_182511.3</u> , <u>BC035789</u> , <u>BC035789.1</u> , <u>BC110843</u> , <u>BX409761</u>
UniProt ID:	<u>Q8IUUK8</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)