

## Product datasheet for **HP219070**

### ZNF714 Human qPCR Primer Pair (NM\_182515)

#### Product data:

|                   |   |
|-------------------|---|
| Product Type:     | qPCR Primer Pairs   |
| Gene ID:          | 148206  |
| Forward Sequence: | GAGATGGTGGATGAATCCCCAG  |
| Reverse Sequence: | ACTCAACCACATTTGCGGAGCC  |
| ACCN:             | <u><a href="#">NM_182515</a></u> , <u><a href="#">NM_182515.1</a></u> , <u><a href="#">NM_182515.2</a></u> , <u><a href="#">NM_182515.3</a></u> , <u><a href="#">BC022527</a></u> , <u><a href="#">BC037782</a></u> , <u><a href="#">BC050468</a></u> ,<br><u><a href="#">BC064985</a></u> , <u><a href="#">BM788748</a></u> , <u><a href="#">NM_182515.4</a></u> |
| UniProt ID:       | <u><a href="#">Q96N38</a></u>   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.                      |
| Storage:          | Store at -20°C.   |
| Stability:        | The primer mix is stable for one year from date of shipping.  |



[View online »](#)