

Product datasheet for **HP219097**

TMED4 Human qPCR Primer Pair (NM_182547)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	222068
Forward Sequence:	TCTCTTCGCTGGTGGCAAACCTG
Reverse Sequence:	GTCCTTCTGAATCTGTTCCACC
ACCN:	<u>BC025705</u> , <u>NM_182547</u> , <u>NM_182547.1</u> , <u>NM_182547.2</u> , <u>NM_182547.3</u> , <u>BC052641</u> , <u>BC052641.1</u> , <u>BC035467</u> , <u>BC053583</u> , <u>BC057851</u> , <u>BM857430</u> , <u>BP872004</u> , <u>BX461603</u> , <u>BX647965</u>
UniProt ID:	<u>Q7Z7H5</u>
Synonyms:	ERS25; GMP25iso; HNLF; p24a3; p24alpha3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)