

## Product datasheet for **HP219099**

### HLA-DQB2 Human qPCR Primer Pair (NM\_182549)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3120
Forward Sequence:	GTGTGCAGACACAACACTACGAGG
Reverse Sequence:	TCACTGAGCAGACCAGCAGGTT
ACCN:	<a href="#">NM_182549</a> , <a href="#">NM_182549.1</a> , <a href="#">BC031995</a>
UniProt ID:	<a href="#">P05538</a>
Synonyms:	HLA-DQB1; HLA-DXB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)