

Product datasheet for **HP219165**

DPY19L2P2 Human qPCR Primer Pair (NM_182634)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	349152
Forward Sequence:	TGGACACCACCTCTCCGTGAAA
Reverse Sequence:	GCAACATTGGAAAGACAGAGTGC
ACCN:	<u>NM_182634</u> , <u>NM_182634.2</u> , <u>BC023575</u> , <u>BC033022</u> , <u>BC045595</u> , <u>BC068601</u> , <u>BC094802</u> , <u>BC150309</u> , <u>BX648212</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)